

QuickTalker Freestyle™

Latitude for QuickTalker Freestyle Mounting System

Item #10000023

Includes:

- Quick release mounting cradle for QuickTalker Freestyle
- 3 tube mounting arm - 30cm, 40cm, and 25cm tubes
- Super Clamp mounting bracket



Mounting Cradle

Easily attach or remove QuickTalker Freestyle

Flexibility

Each joint rotates and allows tubes to move up and down for maximum flexibility

Easy In, Easy Out

Entire mount swivels to allow the user to enter or exit the wheelchair

Super Clamp

Our most versatile clamp for mounting to different surfaces

AbleNet®

quicktalkerfreestyle.com
(651) 294-2211

Latitude for QuickTalker Freestyle Ordering Guide

The following details are required to fulfill your Latitude for QuickTalker Freestyle order

1. Wheelchair Manufacturer:
2. Wheelchair Model #:
3. Please provide two photos showing where on the wheelchair you intend to attach the Latitude QTF mounting arm.



One zoomed out



One zoomed in

4. What is the shape of the tubing where you intend to attach the Latitude QTF mounting arm?

- Circle
- Oval
- Square / Rectangle
- Other

5. Please provide the dimensions of the tubing where you intend to attach the Latitude QTF mounting arm.

For circle or oval tubing, use a soft measuring tape to get the outside circumference of the tubing. If you do not have a soft measuring tape, you can use a piece of string, like a shoe lace to wrap around the tubing, and then measure the length of string with a ruler.

Circle or Oval tubing outside circumference:

Square /Rectangle tubing: Length

Width

PLEASE RETURN COMPLETED FORM TO
fundingservices@ablenetinc.com